

BLOOD PROCESSING FORM I: Heparin Tubes (Orange Top Vacutainers)
Kenya National Micronutrient Survey - 2011

Team # _____

Date _____ / _____ / 2011
Day Month

Cluster # _____

Time samples arrived at cluster lab _____ : _____
Hour Minute

	Target Group (e.g. WRA)	Label ID # (####)	Blood volume at Arrival (mL)	Time of Centrifuge	Time start aliquoting	Haem-olysed <i>Circle</i>	PCR tube <i>circle</i>	PCR ID # (####)	B12 cryovial <i>circle</i>	Zinc cryovial <i>circle</i>	Hbpath cryovial <i>circle</i>	Storage cryovial <i>Circle</i>	Box #		Comments
													PCR tubes	Cryovials	
1.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
2.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
3.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
4.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
5.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
6.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			
7.				____ : ____ Hour Minute	____ : ____ Hour Minute	Y N	Y N		Y N	Y N	Y N	Y N			

Time processed samples were placed in -20°C freezer _____ : _____
Hour Minute

Signature of Cluster Laboratory Technician

ID of Cluster Laboratory Technician